

PDCD1 Human

Programmed Cell Death 1 Human Recombinant

RCP1604

Product Overview

Name	PDCD1 Human
Accession(Primary)	Q15116
Immunogen	Anti-human PDCD1 mAb, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human PDCD1 amino acids 21-167 purified from E. coli.
Description	Programmed Cell Death 1 Human Recombinant

Target information (Q15116)

Gene Name

About Gene Symbol