

PDCD1 Human, Active

Programmed Cell Death 1 Human Recombinant, Active

RCP1605

Product Overview

Name	PDCD1 Human, Active
Accession(Primary)	Q15116
Immunogen	Anti-human PDCD12 mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human PDCD12 amino acids 254-362 purified from E. coli.
Description	Programmed Cell Death 1 Human Recombinant, Active

Target information (Q15116)

Gene Name

About Gene Symbol