

PDCD1LG2 Human, Sf9

Programmed Cell Death 1 Ligand 2 Human Recombinant, Sf9

RCP1607

Product Overview

Name	PDCD1LG2 Human, Sf9
Accession(Primary)	Q9BQ51
Immunogen	Anti-human PDCD4 mAb is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human PDCD4 amino acids 1-469 purified from E. coli.
Description	Programmed Cell Death 1 Ligand 2 Human Recombinant, Sf9

Target information (Q9BQ51)

Gene Name

About Gene Symbol