

PDCL3 Human

Phosducin-Like 3 Human Recombinant

RCP1613

Product Overview

Name	PDCL3 Human
Accession(Primary)	Q9H2J4
Description	Phosducin-Like 3 Human Recombinant

Target information (Q9H2J4)

Gene Name

About Gene Symbol

PDCL3 (Phosducin Like 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.