

PEA15 Human

Phosphoprotein Enriched in Astrocytes 15 Human Recombinant

RCP1619

Product Overview

Name	PEA15 Human
Accession(Primary)	Q15121
Immunogen	Anti-human PEBP1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PEBP1 protein, 1-187 amino acids purified from E. coli.
Description	Phosphoprotein Enriched in Astrocytes 15 Human Recombinant

Target information (Q15121)

Gene Name

About Gene Symbol