

PGLYRP1 Human, HEK

Peptidoglycan Recognition Protein 1 Human Recombinant, HEK

RCP1638

Product Overview

Name	PGLYRP1 Human, HEK
Accession(Primary)	O75594
Immunogen	Anti-human PGM1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PGM1 protein 1-562 amino acids purified from E. coli.
Description	Peptidoglycan Recognition Protein 1 Human Recombinant, HEK

Target information (O75594)

Gene Name

About Gene Symbol