

PLSCR3 Human

Phospholipid Scramblase 3 Human Recombinant

RCP1668

Product Overview

Name	PLSCR3 Human
Accession(Primary)	Q9NRY6
Description	Phospholipid Scramblase 3 Human Recombinant

Target information (Q9NRY6)

Gene Name

About Gene Symbol

PLSCR3 (Phospholipid Scramblase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.