

PSIP1 Human

PC4 And SFRS1 Interacting Protein 1 Human Recombinant

RCP1733

Product Overview

Name	PSIP1 Human
Accession(Primary)	O75475
Immunogen	Anti-human PSMA mAb is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human PSMA amino acids 117-351 purified from E. coli.
Description	PC4 And SFRS1 Interacting Protein 1 Human Recombinant

Target information (O75475)

Gene Name	
About Gene Symbol	PSIP1 (PC4 And SRSF1 Interacting Protein 1) is a transcriptional coactivator involved in adaptive immunity. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.