

RAC1 Human, His

Ras-Related C3 Botulinum Toxin substrate 1 Human Recombinant, His Tag

RCP1782

Product Overview

Name	RAC1 Human, His
Accession(Primary)	P63000
Immunogen	Anti-human RAC2, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human RAC2 amino acids 1-189 purified from E. coli.
Description	Ras-Related C3 Botulinum Toxin substrate 1 Human Recombinant, His Tag

Target information (P63000)

Gene Name

About Gene Symbol