

RANBP1 Human

RAN Binding Protein 1 Human Recombinant

RCP1804

Product Overview

Name	RANBP1 Human
Accession(Primary)	P43487
Immunogen	Anti-human RANBP3 mAb, clone PAT12E11AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human RANBP3 protein (235-445 amino acids) purified from E. coli.
Description	RAN Binding Protein 1 Human Recombinant

Target information (P43487)

Gene Name

About Gene Symbol