

S100A10 Mouse

Recombinant Mouse S100 Calcium Binding Protein A10

RCP1949

Product Overview

Name	S100A10 Mouse
Accession(Primary)	P08207
Immunogen	Anti-human s100A11 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human s100A11 protein 1-105 amino acids purified from E. coli.
Description	Recombinant Mouse S100 Calcium Binding Protein A10

Target information (P08207)

Gene Name

About Gene Symbol