

S100A3 Mouse

S100 Calcium Binding Protein A3 Mouse Recombinant

RCP1958

Product Overview

Name	S100A3 Mouse
Accession(Primary)	P62818
Immunogen	Anti-human s100A4 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human s100A4 protein 1-101 amino acids purified from E. coli.
Description	S100 Calcium Binding Protein A3 Mouse Recombinant

Target information (P62818)

Gene Name

About Gene Symbol