

SAR1A Human

GTP-Binding Protein SAR1A Human Recombinant

RCP1989

Product Overview

Name	SAR1A Human
Accession(Primary)	Q9NR31
Immunogen	Anti-human SAR1B mAb, clone PAT1C7A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human MRRF protein 1-198 amino acids purified from E. coli.
Description	GTP-Binding Protein SAR1A Human Recombinant

Target information (Q9NR31)

Gene Name**About Gene Symbol**