

SERPINB4 Human

Serpin Peptidase Inhibitor, Clade B Member 4 Human Recombinant

RCP2047

Product Overview

Name	SERPINB4 Human
Accession(Primary)	P48594
Immunogen	Anti-human SERPINB5 mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human SERPINB5 amino acids 1-375 purified from E. coli.
Description	Serpin Peptidase Inhibitor, Clade B Member 4 Human Recombinant

Target information (P48594)

Gene Name

About Gene Symbol