

SGCB Human

Sarcoglycan Beta Human Recombinant

RCP2077

Product Overview

Name	SGCB Human
Accession(Primary)	Q16585
Immunogen	Anti-human SGCD mAb, clone PAT19G8AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SGCD protein 57-289 amino acids purified from E. coli.
Description	Sarcoglycan Beta Human Recombinant

Target information (Q16585)

Gene Name**About Gene Symbol**