

SMAD3 Human

Mothers Against Decapentaplegic Homolog 3 Human Recombinant

RCP2130

Product Overview

Name	SMAD3 Human
Accession(Primary)	P84022
Description	Mothers Against Decapentaplegic Homolog 3 Human Recombinant

Target information (P84022)

Gene Name

About Gene Symbol

SMAD3 (SMAD Family Member 3) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.