

SMAD4 Human

Mothers Against Decapentaplegic Homolog 4 Human Recombinant

RCP2131

Product Overview

Name	SMAD4 Human
Accession(Primary)	Q13485
Description	Mothers Against Decapentaplegic Homolog 4 Human Recombinant

Target information (Q13485)

Gene Name

About Gene Symbol

SMAD4 (Mothers Against Decapentaplegic Homolog 4) is a transcription factor involved in TGF- β signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

References for protein:

Title :Transforming Growth Factor β -Induced Reactivation of Epstein-Barr Virus Involves Multiple Smad-Binding Elements Cooperatively Activating Expression of the Latent-Lytic Switch BZLF1 Gene .
Publication :Published ahead of print 18 May 2011, doi: 10.1128/JVI.01197-10 J. Virol. August 2011 vol. 85 no. 15 7836-7848 Link : SMAD4 prospec publication