

SNAP23 Human

Synaptosomal-associated protein 23 Human Recombinant

RCP2137

Product Overview

Name	SNAP23 Human
Accession(Primary)	O00161
Immunogen	Anti-human SNAP25 mAb is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with recombinant human SNAP25 amino acids 1-206 purified from E. coli.
Description	Synaptosomal-associated protein 23 Human Recombinant

Target information (O00161)

Gene Name

About Gene Symbol