

SNCA A53T Human

Alpha Synuclein A53T Human Recombinant

RCP2149

Product Overview

Name	SNCA A53T Human
Accession(Primary)	P37840
Immunogen	Anti-human α -Synuclein mAb is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human α -Synuclein amino acids 61-95 purified from E. coli.
Description	Alpha Synuclein A53T Human Recombinant

Target information (P37840)

Gene Name**About Gene Symbol**