

SNRPG Human

Small Nuclear Ribonucleoprotein Polypeptide G Human Recombinant

RCP2177

Product Overview

Name	SNRPG Human
Accession(Primary)	P62308
Immunogen	Anti-human SNTA1 mAb, clone PAT1E1A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SNTA1 protein 1-505 amino acids purified from E. coli.
Description	Small Nuclear Ribonucleoprotein Polypeptide G Human Recombinant

Target information (P62308)

Gene Name

About Gene Symbol