

SOD Human, 15N

Superoxide Dismutase, 15N Human Recombinant

RCP2188

Product Overview

Name	SOD Human, 15N
Accession(Primary)	P00441
Immunogen	Anti-human SOD1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SOD1 protein 1-154 amino acids purified from E. coli.
Description	Superoxide Dismutase, 15N Human Recombinant

Target information (P00441)

Gene Name**About Gene Symbol**