

## SORBS3 Human

*Sorbin And SH3 Domain Containing 3 Human Recombinant*

RCP2193

### Product Overview

|                           |                                                                                                                                                                                                                           |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>               | SORBS3 Human                                                                                                                                                                                                              |
| <b>Accession(Primary)</b> | O60504                                                                                                                                                                                                                    |
| <b>Immunogen</b>          | Anti-human SORD mAb, clone PAT10F4AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SORD protein 1-357 amino acids purified from E. coli. |
| <b>Description</b>        | Sorbin And SH3 Domain Containing 3 Human Recombinant                                                                                                                                                                      |

### Target information (O60504)

**Gene Name**

**About Gene Symbol**