

SSR2 Human

Signal Sequence Receptor, Beta Human Recombinant

RCP2248

Product Overview

Name	SSR2 Human
Accession(Primary)	P43308
Immunogen	Anti-human SSR4 mAb, clone PAT26G5AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SSR4 protein 24-144 amino acids purified from E. coli.
Description	Signal Sequence Receptor, Beta Human Recombinant

Target information (P43308)

Gene Name

About Gene Symbol