

STMN4 Human

Stathmin Like-4 Human Recombinant

RCP2269

Product Overview

Name	STMN4 Human
Accession(Primary)	Q9H169
Immunogen	Anti-human STOM mAb, clone PAT33F5AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human STOM protein 55-288 amino acids purified from E. coli.
Description	Stathmin Like-4 Human Recombinant

Target information (Q9H169)

Gene Name**About Gene Symbol**