

Streptavidin (37-159), His

Streptavidin (37-159 a.a) Recombinant, His Tag

RCP2274

Product Overview

Name	Streptavidin (37-159), His
Accession(Primary)	P22629
Immunogen	Anti Streptavidin mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant Streptavidin amino acids 25-183 purified from E. coli.
Description	Streptavidin (37-159 a.a) Recombinant, His Tag

Target information (P22629)

Gene Name**About Gene Symbol**