

SUMO2 Human

Small Ubiquitin-Related Modifier 2 Human Recombinant

RCP2295

Product Overview

Name	SUMO2 Human
Accession(Primary)	P61956
Immunogen	Anti-human SUMO2/3 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human SUMO2 amino acids 1-93 purified from E. coli.
Description	Small Ubiquitin-Related Modifier 2 Human Recombinant

Target information (P61956)

Gene Name

About Gene Symbol