

TBCC Human

Tubulin Folding Cofactor C Human Recombinant

RCP2333

Product Overview

Name	TBCC Human
Accession(Primary)	Q15814
Immunogen	Anti-human TBCEL mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TBCEL protein 1-424 amino acids purified from E. coli.
Description	Tubulin Folding Cofactor C Human Recombinant

Target information (Q15814)

Gene Name**About Gene Symbol**