

TCP1 Human

T-Complex 1 Human Recombinant

RCP2348

Product Overview

Name	TCP1 Human
Accession(Primary)	P17987
Immunogen	Anti-human TDG mAb, clone PAT2F7A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TDG protein 1-410 amino acids purified from E. coli.
Description	T-Complex 1 Human Recombinant

Target information (P17987)

Gene Name	
About Gene Symbol	TCP1 (T-Complex 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.