

TFRC Human, SF9

Transferrin Receptor Human Recombinant, Sf9

RCP2359

Product Overview

Name	TFRC Human, SF9
Accession(Primary)	P02786
Description	Transferrin Receptor Human Recombinant, Sf9

Target information (P02786)

Gene Name

About Gene Symbol

TFRC (Transferrin Receptor) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.