

TLR3 Human

Toll Like Receptor 3 Human Recombinant

RCP2388

Product Overview

Name	TLR3 Human
Accession(Primary)	O15455
Immunogen	Anti-human TLR7 mAb, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with recombinant human TLR7 amino acids 451-500 purified from E. coli.
Description	Toll Like Receptor 3 Human Recombinant

Target information (O15455)

Gene Name**About Gene Symbol**