

TMEFF1 Human

TMEFF1 Human Recombinant

RCP2390

Product Overview

Name	TMEFF1 Human
Accession(Primary)	Q8IYR6
Immunogen	Anti-human TMEFF2 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human TMEFF2 amino acids 41-150 purified from E. coli.
Description	TMEFF1 Human Recombinant

Target information (Q8IYR6)

Gene Name

About Gene Symbol