

TNIP1 Human

TNFAIP3 Interacting Protein 1 Human Recombinant

RCP2395

Product Overview

Name	TNIP1 Human
Accession(Primary)	Q15025
Immunogen	Anti-human TNNI1 mAb, clone PAT36E7A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TNNI1 protein.
Description	TNFAIP3 Interacting Protein 1 Human Recombinant

Target information (Q15025)

Gene Name

About Gene Symbol