

TOM1L2 Human

Target of Myb1-Like 2 Human Recombinant

RCP2409

Product Overview

Name	TOM1L2 Human
Accession(Primary)	Q6ZVM7
Immunogen	Anti-human TOMM20 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TOMM20 protein 25-145 amino acids purified from E. coli.
Description	Target of Myb1-Like 2 Human Recombinant

Target information (Q6ZVM7)

Gene Name

About Gene Symbol