

TP53AIP1 Human

Tumor Protein P53 Regulated Apoptosis Inducing Protein 1 Human Recombinant

RCP2413

Product Overview

Name	TP53AIP1 Human
Accession(Primary)	Q9HCN2
Immunogen	Anti-human TP53I3 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TP53I3 amino acids 1-332 purified from E. coli.
Description	Tumor Protein P53 Regulated Apoptosis Inducing Protein 1 Human Recombinant

Target information (Q9HCN2)

Gene Name

About Gene Symbol