

TSFM Human

Ts Translation Elongation Factor Mitochondrial Human Recombinant

RCP2458

Product Overview

Name	TSFM Human
Accession(Primary)	P43897
Description	Ts Translation Elongation Factor Mitochondrial Human Recombinant

Target information (P43897)

Gene Name

About Gene Symbol

TSFM (Ts Translation Elongation Factor, Mitochondrial) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.