

TTC33 Human

Tetratricopeptide Repeat Domain 33 Human Recombinant

RCP2467

Product Overview

Name	TTC33 Human
Accession(Primary)	O96033
Immunogen	Anti-human TTR mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human TTR amino acids 21-147 purified from E. coli.
Description	Tetratricopeptide Repeat Domain 33 Human Recombinant

Target information (O96033)

Gene Name

About Gene Symbol