

TUBB3 Human

Tubulin, Beta 3 Class III Human Recombinant

RCP2468

Product Overview

Name	TUBB3 Human
Accession(Primary)	Q13509
Immunogen	Anti-human TUBG1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TUBG1 protein 1-451 amino acids purified from E. coli.
Description	Tubulin, Beta 3 Class III Human Recombinant

Target information (Q13509)

Gene Name	
About Gene Symbol	TUBB3 (Tubulin Beta 3 Class III) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.