

TXN1, His

Thioredoxin Recombinant, His Tag

RCP2478

Product Overview

Name	TXN1, His
Accession(Primary)	P0AA25
Immunogen	Anti-human TXN2 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TXN2 protein 60-166 amino acids purified from E. coli.
Description	Thioredoxin Recombinant, His Tag

Target information (P0AA25)

Gene Name**About Gene Symbol**