

VAMP2 Mouse

Synaptobrevin-2 Recombinant Mouse

RCP2534

Product Overview

Name	VAMP2 Mouse
Accession(Primary)	P63044
Immunogen	Anti-human VAMP3 mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human VAMP3 amino acids 1-77 purified from E. coli.
Description	Synaptobrevin-2 Recombinant Mouse

Target information (P63044)

Gene Name**About Gene Symbol**