

VASP Human

Vasodilator-Stimulated Phosphoprotein Human Recombinant

RCP2542

Product Overview

Name	VASP Human
Accession(Primary)	P50552
Immunogen	Anti-human Vaspin mAb, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human Vaspin amino acids 21-414 purified from E. coli.
Description	Vasodilator-Stimulated Phosphoprotein Human Recombinant

Target information (P50552)

Gene Name	
About Gene Symbol	VASP (Vasodilator Stimulated Phosphoprotein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.