

Vimentin Human

Vimentin Human Recombinant

RCP2553

Product Overview

Name	Vimentin Human
Accession(Primary)	P08670
Immunogen	Anti-human Vimentin mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Vimentin protein 1-466 amino acids purified from E. coli.
Description	Vimentin Human Recombinant

Target information (P08670)

Gene Name	VIM (Vimentin)
About Gene Symbol	VIM (Vimentin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.