

VSTM2L Human

V-Set and Transmembrane Domain Containing 2 Like Human Recombinant

RCP2566

Product Overview

Name	VSTM2L Human
Accession(Primary)	Q96N03
Immunogen	Anti-human VTA1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human VTA1 protein 1-307 amino acids purified from E. coli.
Description	V-Set and Transmembrane Domain Containing 2 Like Human Recombinant

Target information (Q96N03)

Gene Name

About Gene Symbol