

ANG Human, Sf9

ANG Human Recombinant, Sf9

RCP2624

Product Overview

Name	ANG Human, Sf9
Accession(Primary)	P03950
Immunogen	Anti-human ANGPTL3 mAb, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human ANGPTL3 amino acids 243-460 purified from E. coli.
Description	ANG Human Recombinant, Sf9

Target information (P03950)

Gene Name

About Gene Symbol