

CHGA Human, GST

Chromogranin A Human Recombinant, GST Tag

RCP2637

Product Overview

Name	CHGA Human, GST
Accession(Primary)	P10645
Immunogen	Anti-human CHI3L1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human CHI3L1 amino acids 22-383 purified from E. coli.
Description	Chromogranin A Human Recombinant, GST Tag

Target information (P10645)

Gene Name	
About Gene Symbol	CHGA (Chromogranin A) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.