

S100A8 Human

*S100 Calcium Binding Protein A8 Human Recombinant
RCP1968*

Product Overview

Name S100A8 Human

Description

S100 Calcium Binding Protein A8 Human Recombinant

Accession (Primary) [P05109](#)

Synonyms

Calgranulin A, MRP8, CAGA, CGLA, CFAG, Protein S100-A8, S100 calcium-binding protein A8, Migration inhibitory factor-related protein 8, MRP-8, p8, Cystic fibrosis antigen, Leukocyte L1 complex light chain, Calprotectin L1L subunit, Urinary stone protein band A, S100A8, MIF, NIF, L1Ag, CP-10, MA387, 60B8AG.

Introduction

S100A8 is a part of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a broad range of cells, and participate in the regulation of cellular processes such as cell cycle progression and differentiation. S100A8 plays a role in the inhibition of casein kinase and as a cytokine. S100A8 altered expression is related with cystic fibrosis disease. S100A8 is a calcium-binding protein that has antimicrobial activity against bacteria and fungi. S100A8 is crucial for resistance towards invasion by pathogenic bacteria. S100A8 up-regulates transcription of genes that are under the control of NF-kappa-B. S100A8 plays a role in the development of endotoxic shock in response to bacterial lipopolysaccharide. S100A8 endorses tubulin polymerization and promotes phagocyte migration and infiltration of granulocytes at sites of wounding. S100A8 takes part as a pro-inflammatory mediator in acute and chronic inflammation and up-regulates the release of IL8 and cell-surface expression of ICAM1.

Source

Escherichia Coli.

Formulation

S100A8 was filtered (0.4 µm) and lyophilized in 0.5mg/ml in 20mM Tris and 100mM NaCl, pH 7.5.

Stability

Store lyophilized protein at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles. Reconstituted protein can be stored at 4°C for a limited period of time; it does not show any change after two weeks at 4°C.

Amino acid sequence

MKHHHHHHAS MLTELEKALN SIIDVYHKYS LIKGNFHAVY RDDLKLEET ECPQYIRKKG ADVWFKELDI
NTDGAVNFQEMLTELEKALN SIIDVYHKYS LIKGNFHAVY RDDLKLEET ECPQYIRKKG ADVWFKELDI
NTDGAVNFQE FLILVIKMGV AAHKKSHEES HKE.

Solubility

It is recommended to add deionized water to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

Precautions

S100A8 Human is for research use only and not for use in diagnostic or therapeutic procedures.

Target Information: ([P05109](#))