

## S6K1 (3A2) Mouse mAb

ABZ4085


### Product Overview

|                             |                      |
|-----------------------------|----------------------|
| <b>Name</b>                 | S6K1 (3A2) Mouse mAb |
| <b>Accession(Primary)</b>   | P23443               |
| <b>Host Species</b>         | Mouse                |
| <b>Clonality</b>            | Monoclonal           |
| <b>Isotype</b>              | IgG1                 |
| <b>Species Reactivity</b>   | Human, Mouse, Rat    |
| <b>Recommended Dilution</b> | IHC-1:50-1:100       |


### Target information (P23443)

**Gene Name** RPS6KB1

#### About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using S6K1 (3A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using S6K1 (3A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.