

SAE1 Rabbit mAb

ABZ1782


Product Overview

Name	SAE1 Rabbit mAb
Accession(Primary)	Q9UBE0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100


Target information (Q9UBE0)

Gene Name SAE1

About Gene Symbol



Immunocytochemistry analysis of SAE1 (green) in Hela using SAE1 antibody, and DAPI(blue).



Western blot analysis of SAE1 in HeLa, Jurkat
lysates using SAE1 antibody.