

## SCAI Rabbit mAb


ABZ1992

### Product Overview


|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Name</b>                 | SCAI Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | Q8N9R8                                                                  |
| <b>Host Species</b>         | Rabbit                                                                  |
| <b>Clonality</b>            | Monoclonal                                                              |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                       |
| <b>Calculated MW</b>        | Calculated MW: 70 kDa; Observed MW: 60 kDa                              |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:20-1:50;<br>ICC/IF-1:200-1:500; FC-1:100-1:200 |

### Target information (Q8N9R8)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | SCAI |
| <b>About Gene Symbol</b> |      |



Western blot analysis of SCAI in C6, 3T3 lysates using SCAI antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SCAI antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.