

## SEC61A Rabbit mAb


ABZ1830

### Product Overview


|                             |                                                                               |
|-----------------------------|-------------------------------------------------------------------------------|
| <b>Name</b>                 | SEC61A Rabbit mAb                                                             |
| <b>Accession(Primary)</b>   | Q9H9S3                                                                        |
| <b>Host Species</b>         | Rabbit                                                                        |
| <b>Clonality</b>            | Monoclonal                                                                    |
| <b>Isotype</b>              | IgG                                                                           |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                             |
| <b>Calculated MW</b>        | Calculated MW: 52 kDa; Observed MW: 40 kDa                                    |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200;<br>IP-1:20; ICC/IF-1:50-1:100 |

### Target information (Q9H9S3)


|                          |         |
|--------------------------|---------|
| <b>Gene Name</b>         | SEC61A2 |
| <b>About Gene Symbol</b> |         |



Western blot analysis of SEC61A in K562, C6, 3T3 lysates using SEC61A antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SEC61A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of SEC61A(green)  
in Hela using SEC61A antibody, and DAPI(blue)