


SEH1L Rabbit mAb

ABZ1984

Product Overview

Name	SEH1L Rabbit mAb
Accession(Primary)	Q8R2U0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q8R2U0)

Gene Name Seh1L**About Gene Symbol**

Western blot analysis of SEH1L in K562, C6, 3T3 lysates using SEH1L antibody.

